


















Symbol	Standard	Ref #	Title	Definition
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.1.1	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.1.2	Authorized Representative in the European Community	Indicates the Authorized Representative in the European Community.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.1.3	Date of Manufacture	Symbol for date of manufacture. This symbol is accompanied by a date.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.1.4	Use-by Date	Indicates the date after which the medical device is not to be used.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.1.5	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.1.6	Catalog Number	Indicates the manufacturer's catalog number so that the medical device can be identified.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.2.3	Sterilized Using Ethylene Oxide	Indicates a medical device that has been sterilized using ethylene oxide.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.2.7	Non-Sterile	Indicates a medical device that has not been subjected to a sterilization process.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.2.8	Do Not Use if Package is Damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.4.2	Do Not Reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.4.3	Consult Instructions For Use	Indicates the need for the user to consult the instructions for use.
	ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	5.4.4	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

	21 CFR 801.15(c)(1)(i)F		Prescription only	Indicates a medical device that requires prescription in the United States.
	765/2008/EC 768/2008/EC MDD 93/42/EEC (Articles 4,11,12,17, Annex II)		CE marking	Signifies European technical conformity.
	IEC 60417-5333		Applied Part - Type BF	Indicates a device that comes into contact with the body, having medium or long term contact with the patient.
	BS EN 50419		WEEE (Waste Electrical and Electronic Equipment)	Indicates the device is WEEE per European Community Directive 2012/19/EU
	ISO 7010	M002	Refer to accompanying documents	To signify that the instruction manual/booklet must be read